

BENUTZERINFORMATIONEN/INFORMATION FOR USER BREASEY Med

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


*)

Medical face masks

Medical face masks are mainly used for the protection of others, as they do not close off the face completely. The masks reduce the leakage of liquid droplets from the wearer’s mouth and nose and thus protect others. They are used, for example, to prevent droplets from the patient’s breath from entering open wounds.

Types of medical face mask

Medical face masks are divided into two types according to their bacterial filtration capacity: Type I and Type II, whereby Type II is further differentiated according to whether the mask is also liquid-resistant or not. Table 1 shows the normative performance requirements for the different mask types.

Type I medical face masks should only be used on patients and others to reduce the risk of spreading infection. Type I masks are not intended for use by healthcare professionals in operating rooms or other facilities with similar requirements.

Type II medical face masks intended for use in operating theatres and healthcare facilities with similar requirements are designed to protect the entire work environment. Type II masks are generally approved for use by medical personnel in operating theatres or other medical facilities with similar requirements. These masks mainly protect the patient.

Masks **Type II with the suffix “R”** (Type IIR) are resistant to liquids. They provide spray resistance in accordance with EN 14683:2019+AC:2019, which provides the wearer with protection against saliva, blood and other bodily fluids of foreign persons.

Important information

Medical face masks do not replace respiratory protection. This is not personal protective equipment according to the PPE regulation ((EU) 2016/425). The mask does not protect the wearer from infection. The user must be familiar with the use and handling of the product.

ATTENTION: Medical face masks are for single use and must not be reused. Pathogens attach themselves to the mask and, if reused, could contribute to the spread of germs. Moisturised masks should be changed as they lose most of their protective function.

Do not move or touch the medical face mask while wearing it. It must not hang around the neck of the user at any time.

All serious incidents relating to this product must be notified to the manufacturer and the competent authority of the Member State where the user and/or patient is established.

Storage

The mask must be stored dry in its packaging at room temperature and protected from sunlight. The product must also be protected from extreme temperatures, weather influences and magnetic radiation. If the packaging is damaged, the mask must not be used.

Even if stored properly, the mask can only be used until the expiry date, maximum 5 years after the production date. (a.)

The full declaration of conformity can be accessed via the following link: www.skylotec.com/downloads

Properties of the mask

- Non-sterile medical face mask (b.)
- Not reusable (c.)
- It is a medical device (d.)
- Identify the outside of the mask. The stamping “out” stands for “outside” and is located at the bottom edge on the left half of the face. (e.)
- Tested according to EN 14683:2019+AC:2019

Performance requirements for medical face masks (Table 1)

TYP	Bacterial Filter performance	Respiratory Resistance	Pressure of the spray resistance	Microbiological Clarity
TYP I	> 95 %	< 40 Pa/cm ²	-	< 30 KBE/g
TYP II	> 98 %	< 40 Pa/cm ²	-	< 30 KBE/g
TYP IIR	> 98 %	< 60 Pa/cm ²	> 16,0 kPa	< 30 KBE/g

*) Product/original can deviate from the illustration
Subject to technical changes

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© SKYLOTEC