



BREASEY Made in Germany

BREASEY MNS IIR



DE Atmungsaktiv & geruchsfrei
EN Breathable and odorless
FR Respirable et inodore
ES Transpirable e inodoro
IT Traspirante e inodore



Weiches Material
Soft material
Matériau souple
Material suave
Materiale morbido



3-Lagen-Technologie
3 layers technology
Technologie à 3 couches
Tecnologia de 3 capas
Tecnologia a 3 strati



Anpassbarer Nasenbügel
Adjustable nose clip
Arête nasale réglable
Punte nasale regolabile
Ponte nasale regolabile



ElasticEarband

DE Hochelastisches Ohrband
EN Highly elastic earloop
FR Bandeau d'oreille très élastique
ES Banda auricular muy elástica
IT Banda dell'orecchio altamente elastica

DE	EN	FR	ES	IT
Norm	Standard	Standard	Norma	Norma
EN 14683:2019+AC:2019				
Filterleistung	Filter performance	Performance des filtres	Rendimiento del filtro	Prestazioni del filtro
≥ 99 %				
Größe	Size	Taille	Talla	Dimensione
175 mm x 97 mm				
Standard-Farbe*	Standard colour*	Couleur standard*	Color estándar*	Colore standard*
medizinisches Blau	medical blue	bleu médical	azul médico	blu medico
Lagen	Layer	Sites	Ubicaciones	Sedi
3				
Nasenbügel	Nose clip	Pince-nez	Pinza para la nariz	Clip per il naso
✓				
Ohrband	Earloop	Bandeau auriculaire	Banda para el oído	Auricolare
ElasticEarband				
Kopfband	Headband	Bandeau de tête	Cinta para la cabeza	Archetto
-				
Mindesthaltbarkeit /Jahre	Minimum shelf life / years	Durée minimale de conservation / ans	Vida útil mínima / años	Durata minima di conservazione / anni
5				
Material made in Germany	Material made in Germany	Matériel fabriqué en Allemagne	Material fabricado en Alemania	Materiale prodotto in Germania
PP non-woven, Meltblown	PP non-woven, Meltblown	PP filé ; non-tissé soufflé à Chaud	PP meltblown no tejido	PP non-tessuto, Polimero fuso
Masken einzeln verpackt	Masks individually packed	Masques individuels emballés	Máscaras individuales embaladas	Maschere individualmente imballate
-				
Individualisierung	Individualisation	Individualisation	Individualización	Individualizzazione
Auf Anfrage möglich	Possible on request	Possible sur demande	Posible bajo petición	Possibile su richiesta

DE	EN	FR	ES	IT
Artikelnummer	Article number	Numéro d'article	Número de artículo	Numero di articolo
MA-104-710-50				
Masken pro Box	Masks per box	Masques par boîte	Máscaras por caja	Maschere per scatola
50				
Box Abmessung	Box size	Taille de la boîte	Tamaño de la caja	Dimensione della scatola
192 x 104 x 131 mm				
Box Gewicht	Box weight	Poids de la boîte	Peso de la caja	Peso della scatola
0,21 kg				
Box GTIN	Box GTIN	Boîte GTIN	Caja GTIN	Box GTIN
4030281010845				
Boxen pro Umkarton	Boxes per carton	Boîtes par carton externe	Cajas por cartón exterior	Scatole per cartone esterno
15				
Masken pro Umkarton	Masks per carton	Masques par carton	Máscaras por caja	Maschere per cartone
750				
Umkarton Gewicht	Carton weight	Poids du carton	Peso de la caja	Peso del cartone
3,095 kg				
Masken pro Palette	Masks per palette	Masques par palette	Máscaras por paleta	Maschere per palette
18000				
Umkartons pro Palette	Cartons per pallet	Cartons par palette	Cajas por palet	Cartoni per pallet
24				
Palette Gewicht	Pallet weight	Poids des palettes	Peso de la paleta	Peso del pallet
130 kg				
Palette Abmessung	Pallet dimension	Dimension de la palette	Dimensión de la paleta	Dimensione del pallet
1200 x 800 x 1849 mm				
Zolltarifnummer	Tariff number	Numéro de tarif	Número de tarifa	Numero di tariffa
90200000				

* weitere Farben auf Anfrage/ Other colours on request/ Autres couleurs sur demande/ Otros colores a petición/ Altri colori su richiesta